

Customer Name _____ Customer No. _____ Purchase Order No. _____ Fax _____
 Billing Address _____ Patient Name _____
 Shipping Address _____ Order Date _____ Measured By _____
 Telephone _____ Email _____ Shipping Method _____
 Credit Card Info _____

upper leg kit: Check the box next to the corresponding length and width.

*cG cannot be more than 30cm larger than cE2.

	Length			Width		Quantity
	Short	Standard	Long	Regular	Wide	
Left	ℓ D1-G < 37cm	ℓ D1-G 37 - 41cm	ℓ D1-G > 41	cG < 80cm cE2 < 65cm	cG < 125cm cE2 < 110cm	
Right						

Each upper leg kit includes: • 1 reduction upper leg component • 1 Built-In-Tension system guide card
 • 1 reduction shelf strap • 1 paper measuring tape
 • 6 hook and loop stays • 1 direction for use.

knee kit: (the reduction knee kit is only available in one length and width)

*cE2 cannot be more than 30cm larger than cC.

	Length	Width	Quantity
Left	30cm	cE2 < 110cm cE < 102cm cC < 80cm	
Right			

Each knee kit includes: • 1 reduction knee component • 1 paper measuring tape
 • 2 hook and loop knee spines • 1 direction for use
 • 1 reduction shelf strap

lower leg kit: Check the box next to the corresponding length and width.

*cC cannot be more than 30cm larger than cB.

	Length			Width		Quantity
	Short	Standard	Long	Regular	Wide	
Left	ℓ mall-D1 < 32cm	ℓ mall-D1 32 - 36cm	ℓ mall-D1 > 36cm	cC < 65cm cB < 55cm	cC < 80cm cB < 65cm	
Right						

Each lower leg kit includes: • 1 reduction lower leg component • 1 Built-In-Tension system guide card
 • 1 pair reduction kit undersleeve leg • 1 paper measuring tape
 • 1 reduction shelf strap • 1 direction for use
 • 6 hook and loop stays

foot options: Check the box next to the corresponding length and width.

*foot options are not included with the kit and must be ordered separately.


	compression ankle		pac band™		Single band AFW®	Quantity
	Standard	Large	X-Large			
Left	cA1 19-36 cm	cA1 22-41 cm	cA1 < 44cm		cA1 > 41cm	
Right	cB 19-27 cm	cB 25-50 cm	cB < 56cm			

compression ankle available in Standard and Large. pac band™ available in Standard, Large, and X-Large.

additional ancillary items:

reduction kit undersleeve leg Quantity (Pairs)	reduction kit shelf strap Quantity (Each)

G	Groin
E2	15cm above E
E1	5cm above E
E	Center of patella
D1	Knee crease
D	Slightly below knee
C	Widest part of calf
B	Narrowest part of ankle above malleolus
A1	Middle of foot
A	Ball of foot
mall	Center of malleolus
ℓ	length in a straight line
c	circumference




Customer Name _____ Customer No. _____ Purchase Order No. _____ Fax _____
 Billing Address _____ Patient Name _____
 Shipping Address _____ Order Date _____ Measured By _____
 Telephone _____ Email _____ Shipping Method _____
 Credit Card Info _____

arm kit: Check the box next to the corresponding length and width.

*** If cB is > 25cm hand wrap will not fit, though reduction arm component may still be used.**

	Length		Width		Quantity
	Standard	Long	Regular	Wide	
Left	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Right	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	



Each arm kit includes:

- 1 reduction arm component
- 1 pair undersleeve arm
- 1 customizable hand wrap
- 1 Built-In-Tension guide card
- 1 direction for use

Arm Measurement Guide	
G	Axilla
F	Midpoint of E and G
E	Elbow crease
D	Midpoint of C and E
C	Wrist
B	Base of thumb webbing
	ℓ = length in a straight line
	c = circumference

additional ancillary items:

undersleeve arm standard Standard <50cm max circ. Quantity (Each)	undersleeve arm wide Wide <80cm max circ. Quantity (Each)	reduction kit shelf strap Quantity (Each)
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>